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(Human Placenta
Growth Factor)

DRG

PLGF
Freies Estriol
HCG
Freies beta-HCG
AFP
PAPP-A
HPL

DRG

DRG ELISA

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A New Marker for Preeclampsia

In 5-8 percent of all pregnancies, Preeclampsia causes an increase in the mother's blood pressure, excretion of protein, and edemas – in the majority of cases, this occurs in the last trimester of pregnancy.

These clinical conditions, together termed Preeclampsia, are a risk for mother and child which can make it necessary to induce premature birth.

Studies about Preeclampsia have shown that important growth factors are affected. This causes malnutrition of placenta and fetus.

It is now possible to measure the concentration of one of the growth factors involved – PLGF – in blood.

*The PLGF ELISA allows the early detection (**15th week of pregnancy**) of the developmental risk factors of preeclampsia and is therefore an important new tool in the health protection of mother and child.*

PLGF Elisa EIA-4529

Inkubationszeit/Incubation time: 2,5 h
Standardbereich/Standard range: 0-1000 pg/mL
Probe/Sample: 25 µL Serum
Empfindlichkeit/Sensitivity: 1,06 pg/mL
Reagenzien sind gebrauchsfertig/Ready to use reagents
Kitinterne Kontrollen/Internal kit controls

Assay characteristics

CV Intra assay: 1,68 % - 2,83 %
CV Inter assay: 4,10 % - 7,00 %
Recovery: 87,0 % - 105,5 %
Linearity: 88,1 % - 112,6 %

Angiogenic factors in normal and preeclamptic pregnancies

Schmidt M¹, Hoffmann B¹, Callies R¹, Janetzko A², Kimmig R¹, Kasimir-Bauer S¹

¹Department of Gynecology and Obstetrics, University of Essen, Germany and ²DRG Instruments GmbH, Marburg, Germany

Background:

Angiogenesis and vascular transformation are important processes in the normal development of the placenta. Abnormal angiogenesis and vascular transformation are considered to be one of the main reasons for preeclamptic pregnancies and intrauterine growth retardation. Placental growth factor (PLGF), a member of the VEGF family, is produced chiefly by the placenta and is a potent angiogenic factor. The availability of a test to predict preeclampsia would be a powerful tool in preventing preeclampsia-induced risks for the mother and fetus. Now, a new PLGF-ELISA kit for routine diagnostics is available, which offers new possibilities in the prediction of preeclampsia.

Materials and methods:

Using solid phase enzyme-linked immunosorbent assay (ELISA) longitudinal serum concentrations of PLGF (DRG Instruments, Marburg, Germany) were measured in normal pregnancies. The values were compared to those of preeclamptic pregnancies.

Results:

The PLGF concentration in the normal pregnancies showed a steady increase with a peak at 28 to 32 weeks and a consistent decline thereafter. The preeclamptic pregnancies showed a significant lower serum concentration of PLGF compared to the non-preeclamptic pregnancies ($p < 0,05$). There was just a slight peak later in pregnancy.


Figure 1: PLGF expression in normal pregnancies in the course of pregnancy


Figure 2: PLGF expression in normal and preeclamptic pregnancies (> 20 weeks)


Figure 3: PLGF at 15 to 18 weeks (Schmidt M et al, Gyn Geb Rundsch 2008, in press)


Figure 4: PLGF in preeclamptic pregnancies


Figure 6: PLGF and probability of preeclampsia


Figure 5: PLGF Cutoff 41,84 pg/mL
Relative risk: 18,4 (95 % CI: 2,40-140,83)
Odds ratio: 30 (95 % CI: 3,21-280,31)

DRG PLGF-ELISA

- All the measured assay characteristics fulfilled the required specifications
- The assay can be considered reliable
- Comfortable using
 - Kit supplies „ready to use reagents“
 - Short incubation period (2.5h)
 - Attractive price

Results

- PLGF expression showed a steady increase beginning at 15 weeks
- PLGF expression is significantly lower in preeclamptic pregnancies
- Aspirin (ASS) has a protective effect starting < 18 weeks

PLGF as a useful marker?

Yes: PLGF expression in early pregnancy is lower in pregnancies with developing preeclampsia

- Helps to identify high risk patients:
- 15-18 weeks: < 42 pg/mL
- Aspirin therapy in early stage of pregnancy could be a strategy to reduce the risk of preeclampsia

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DRG Diagnostics

DRG Instruments GmbH, founded in 1973 by Dr. Geacintov, subsidiary of DRG Intl. Inc., USA, is a diagnostics manufacturer of ELISAS. The DRG Group operates through a network of DRG subsidiaries in Germany, Poland, Russia and the Czech Republic and through distributors worldwide.

ELISAS that perform


DRG develops and manufactures diagnostic ELISA test kits for use in clinical and research laboratories. The experience of our production and management team guarantees to provide high quality products, competitive prices and excellent customer service.

DRG works to DIN EN ISO 9001:2008, ISO 13485:2007 and ISO 13485:2003 under CMDCAS standard, certified by TÜV Rheinland Product Safety GmbH, an indication of our commitment to customer service, quality control and improved health care.

DRG Instruments GmbH, Germany
Frauenbergstraße 18
D-35039 Marburg
Tel. +49 (0) 64 21/170 00,
Fax +49 (0) 64 21/1700 50
Internet: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by

DRG International Inc. USA
1167 U.S. Highway 22 East
Mountainside, N.J. 07092 USA
Phone: +1 (908) 233-2079
Fax +1 (908) 233-0758
Internet: www.drg-international.com
E-mail: corp@drg-international.com