# SALIVA ELISAS

## Speichel Analytik

# DRG

Testosterone  
Progesterone  
Estradiol  
17alpha-OH Progesterone  
DHEA  
Estriol  
Cortisol  
Lactate

# DRG

# DRG ELISAS

## SALIVA Immunoassays Speichel Analytik

**Sensitive ELISAS**  
for detection of free steroids  
and other parameters in saliva  
samples

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### Benefits:

- Non invasive sampling, can be done also by patients
- Measurement of the free (unbound) fraction of steroids
- Much higher diagnostic relevance with reproducible and reliable results, compared to serum analytics
- Simple and patient-friendly measurement of hormone profiles (Diurnal-, monthly **profiles**)
- Easy sampling in babies and kids
- The measurement of steroids in **veterinary** saliva samples is possible without any further pre-treatment

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- Cat. No.: SLV-4188
- Sensitivity: 0.4 pg/mL
  - Standard range: 1 - 100 pg/mL

## Estrone

- Cat. No.: SLV-5228
- High Sensitivity: 4.7 pg/mL
  - Standard range: 10-300 pg/mL

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### Sampling device

As saliva contains only tiny amounts of protein non-polar analytes such as steroids tend to adsorb to various plastic surfaces. Glass would tend to be the best material for saliva sampling and for the most part Polyethylene, although it is still in use in many laboratories worldwide, should be avoided. However there is an inexpensive polypropylene device available from DRG (Sali-Tubes). It does not show any interference in saliva even in Progesterone testing, which is the most critical steroid.

### Blood in saliva

Traces of blood might interfere with the results of saliva testing. It is therefore important that the salivary samples should not show any red colour. In cases of visible blood traces in saliva the sample should be discarded. After rinsing the device several times with tap water a new saliva sample can be taken but only after waiting for approximately 15 minutes.

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
## DRG Diagnostics

DRG Instruments GmbH, founded in 1973 by Dr. Geacintov, subsidiary of DRG Intl. Inc., USA, is a manufacturer of ELISAS. The DRG Group operates through a network of DRG subsidiaries in Germany, Poland, Russia and the Czech Republic and through distributors worldwide.

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## ELISAS that perform


DRG develops and manufactures diagnostic ELISA test kits for use in clinical and research laboratories. The experience of our production and management team guarantees to provide high quality products, competitive prices and excellent customer service.

DRG works to DIN EN ISO 9001:2008, ISO 13485:20012 and ISO 13485:2003 under CMDCAS standard, certified by TÜV Rheinland Product Safety GmbH, an indication of our commitment to customer service, quality control and improved health care.

DRG Instruments GmbH, Germany  
 Frauenbergstraße 18  
 D-35039 Marburg  
 Tel. +49 (0) 64 21/170 00,  
 Fax +49 (0) 64 21/17 00 50  
 Internet: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
 E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by

**DRG International Inc. USA**  
 841 Mountain Avenue  
 Springfield, N.J. 07081 U.S.A.  
 Tel: +1 (9 73) 564-75 55  
 Fax +1 (9 73) 564-75 56  
 E-Mail: [corp@drg-international.com](mailto:corp@drg-international.com)  
[www.drg-international.com](http://www.drg-international.com)