

● CHROMagar™ Candida


For isolation and differentiation
of major clinical-significant *Candida* species

Plate Reading

● *Candida albicans*
→ green

● *Candida tropicalis*
→ metallic blue

● *Candida krusei*
→ pink, fuzzy

For isolation and differentiation of major clinical-significant *Candida* species

Background

Perhaps as a result of progress in medical therapy with a related decrease in patient resistance and increase in the number of susceptible hosts, the incidence of bloodstream *Candida* infections rose through the 1990s but has stabilised in recent years. More commonly, *Candida* species are involved in superficial oropharyngeal and urogenital infections particularly for immuno-depressed people such as the elderly and AIDS victims. "Early diagnosis is essential for early effective management of the patients." WHO Guidelines on Standard Operating Procedures for Laboratory Diagnosis of HIV-Opportunistic Infections)

Although *C.albicans* remains the major species involved, other types such as *C.tropicalis*, *C.krusei* or *C.glabrata* have increased proportionately as new antifungal agents have worked very effectively against *C.albicans*. This shows the importance of an accurate detection for a proper antifungal therapy choice.

Medium Performance

- 1 **HIGH CONTRAST**
Intense colony colouration helping to differentiate species.
- 2 **FAST**
36-48h incubation.
- 3 **PANORAMIC VIEW**
Panoramic view of the *Candida* population mix with the ability to detect major species but also the presence of minor populations.
- 4 **POWERFUL TOOL**
It displays high specificity and sensitivity for 3 of the major *Candida* species at a time.
Specificity and sensitivity:
Candida albicans → Exceed 99%*
Candida tropicalis → Exceed 99%*
Candida krusei → Exceed 99%*
- 5 **LESS WORKLOAD**
With traditional media such as Sabouraud agar, identification of mixed cultures is difficult, whereas CHROMagar™ Candida, just by the colour of the colony, will instantly differentiate various *Candida* species.
- 6 **CONVENIENT**
Facilitates the choice of appropriate and optimised antifungal therapies.

Medium Description

Powder Base	
Total	47.7 g/L
Agar	15.0
Peptone	10.2
Chromogenic mix	22.0
Chloramphenicol	0.5
Storage at 15/30°C - pH: 6.1 ± 0.2.	
Shelf Life	3 years

Usual Samples	skin, sputum, urine, vaginal specimens
Procedure	Direct Streaking. Incubation 36-48h, 30-37°C. Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Quality Control Strains

C. krusei ATCC® 14243 pink and fuzzy
C. albicans ATCC® 60193 green
C. tropicalis ATCC® 1369 metallic blue
C. glabrata ATCC® 2001 mauve
E. coli ATCC® 25922 inhibited

ATCC® is a registered trademark of the American Type Culture Collection


Order References

Please use these product references when contacting your local distributor:

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