CHROMagar™ **Malassezia**

Instructions For Use

Available in several languages

NT-EXT-029

Version 3

ENGLISH

English Version

FRANCAIS

Version Française

ESPAÑOL


Version Español

DEUTSCH

Deutsch Version

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MEDIUM PURPOSE

Chromogenic medium for detection of *Malassezia* spp.

Malassezia is a fungi naturally found on the skin of many animals, and also on humans. On normal healthy skin it does not cause infections, but when the environment of the skin is altered, *Malassezia* species are able to cause several cutaneous diseases as severe dermatitis or otitis (inflammation of the skin or ears respectively).

COMPOSITION

The product is composed of a single powder medium.

Product	=	Pack
Total g/L		56.3g/L
Composition g/L		Agar 15.0 Peptones and extracts 38.0 Chloramphenicol 0.5 Chromogenic mix 2.8
Aspect		Powder Form
STORAGE		15/30°C
FINAL MEDIA pH		5.8 +/- 0.2

CAUTION :
formula to be supplemented with Glycerol at 2g/L and Tween 40 at 10g/L (not provided in the kit).

PREPARATION (Calculation for 1L)

Step 1

Preparation

- Disperse slowly 56.3 g of powder in 1L of purified water.
 - Add 2g of Glycerol and 10g of Tween 40.
 - Stir until agar is well thickened.
 - Heat and bring to boil (100°C) while swirling or stirring regularly.
- DO NOT HEAT TO MORE THAN 100°C. DO NOT AUTOCLAVE AT 121°C.

Warning 1: If using an autoclave, do so without pressure.

Advice 1: For the 100°C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

Step 2

Pouring

- Cool in a water bath at 45-50°C, swirling or stirring gently.
- Pour medium into sterile Petri dishes.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to one month under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.

INOCULATION

Related samples can be processed by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate in aerobic conditions at 30-37°C for 72 hours.

Typical Samples

Clinical and veterinary samples (skin, ear canals etc.)


Direct streaking or spreading technique

CHROMagar™ Malassezia

INTERPRETATION

Microorganism	Typical colony appearance
<i>Malassezia furfur</i>	→ large, pale pink and wrinkled
Other <i>Malassezia</i> species (including <i>M. globosa</i> & <i>M. restricta</i>)	→ mostly pink to purple
<i>C. albicans</i>	→ green
<i>C. tropicalis</i>	→ metallic blue + mauve halo
<i>C. krusei</i>	→ pink, fuzzy
Other yeast species	→ white to mauve

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Specificity and sensitivity for *Malassezia furfur* are close to 100% (Kaneko *et al.* 2007).
- Definite identification requires additional testing.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>M. furfur</i> ATCC® 14521	→ mauve
<i>C. glabrata</i> ATCC® 2001	→ mauve
<i>C. tropicalis</i> ATCC® 1369	→ metallic blue + mauve halo
<i>C. albicans</i> ATCC® 10231	→ pale green
<i>E. coli</i> ATCC® 25922	→ inhibited
<i>C. krusei</i> ATCC® 14243	→ pink, fuzzy

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity

Pack Size	Ordering References	Weight
5000 ml 250 Tests of 20ml	MZ282	Weight: 56.3gr
25 L 1250 Tests of 20ml	MZ283-25	Weight: 281.5gr

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection
NT-EXT-029 V3 / 05-Mars-15

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