E.sakazakii Agar

Instructions For Use
Available in several languages

NT-EXT-070

Version 2

ENGLISH

English Version

FRANCAIS

Version Française

ESPAÑOL

Version Español


DEUTSCH

Deutsch Version

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E.sakazakii agar plate

E.sakazakii Agar

MEDIUM PURPOSE

Chromogenic medium for the detection of *E.sakazakii* (*Cronobacter* spp) in milk or milk products according to ISO/TS 22964 standard requirements.

« *Enterobacter sakazakii* is a gram-negative, non-spore-forming bacterium belonging to the Enterobacteriaceae family <.....> It has been implicated in outbreaks causing meningitis or enteritis, especially in infants. In the few outbreaks reported 20% to >50% of the infants who contracted the disease died. For survivors, severe lasting complications can result including neurological disorders. The outcome related to adult disease seems to be significantly milder.”

Extracts From WHO Microbial Risk Assessment Feb-2007

COMPOSITION

The product is composed of a single powder medium.

Product	=	Pack
Total g/L		30.7 g/L
Composition g/L		Agar 15.0 Casein Peptone 7.0 Yeast extract 3.0 NaCl 5.0 Sodium deoxycholate 0.6 X-α-Glucoside 0.15 Crystal violet 0.002
Aspect		Powder Form
STORAGE		15/30°C
FINAL MEDIA pH		7.0 +/- 0.2

PREPARATION (Calculation for 1L)

Step 1

Preparation

- Disperse slowly 30.7 g of powder base in 1L of purified water.
- Stir until agar is well thickened.
- Heat to 121°C +/- 1°C during 15 min.

Step 2

Pouring

- Cool in a water bath to 44-47°C, swirling or stirring gently.
- Pour medium into sterile Petri dishes.
- Let it solidify and dry. **Appearance of the agar: translucent, purplish.**

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to two months under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.

INOCULATION

Related samples can be processed by direct streaking on the plate, as well as prior appropriate pre-enrichment and enrichment step according ISO/TS 22964 standards.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- From an appropriate enrichment broth, streak a loopful (10µl) onto plate.
- Incubate at 44°C ± 1°C for 24h ± 2h, in aerobic conditions.

Typical Samples

e.g. Milk powder, powdered infant formula, related environmental samples.


Possible enrichment step
Direct streaking
or spreading technique

E.sakazakii Agar

INTERPRETATION

Microorganism	Typical colony appearance
<i>E.sakazakii</i> (<i>Cronobacter</i> spp)	→ green to blue
Other gram negative bacteria	→ inhibited, colourless, light purple and transparent
Gram positive bacteria	→ inhibited

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Some rare strains of *E.sakazakii* may not support a growth at 44°C.
- Definite identification as *E.sakazakii* requires additional testing, for instance the presence of a yellow pigmentation and the biochemical confirmation tests as described in the ISO/TS 22964 Technical Specification.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>E.sakazakii</i> ATCC® 51329	→ green-blue
<i>E.coli</i> ATCC® 25922	→ grey with light blue
<i>E.cloacae</i> ATCC® 35030	→ grey with light blue
<i>E.faecalis</i> ATCC® 29212	→ inhibited
<i>S.aureus</i> ATCC® 25923	→ inhibited
<i>K.pneumoniae</i> ATCC® 13883	→ partially inhibited

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity

 Pack Size

5000 ml

250 Tests
of 20ml

=

Ordering References

CS812

Weight: 153.5gr

Need some
Technical Documents?

Available
for download on
www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection
NT-EXT-070 V2 / 21-Oct-13

CHROMagar
The Chromogenic Media Pioneer

 CHROMagar 4 place du 18 juin 1940
75006 Paris - France
Email: CHROMagar@CHROMagar.com
Tel +33 (0)1.45.48.05.05. Website: www.CHROMagar.com

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