

## ● CHROMagar™ COL-APSE


Click below:



● CHROMagar™  
**COL-APSE**

**For detection of Colistin resistant  
Gram-negative bacteria**

## Plate Reading

-  • COL-R *E. coli*  
→ Dark pink to reddish
-  • COL-R *Klebsiella*, *Enterobacter*, *Citrobacter*  
→ Metallic blue
-  • COL-R *Pseudomonas*  
→ Translucent cream to green
-  • COL-R *Acinetobacter*  
→ Cream, opaque

## For detection of Colistin resistant Gram-negative bacteria

### Background

Polymyxin E (colistin) and B are increasingly used as antimicrobials in the treatment of multi-drug resistant bacterial infections. Polymyxin resistance, although intrinsic in Gram-positive and some Gram-negative species (*Proteus*, *Morganella*, *Serratia*), is now a problem in a number of other pathogens (*Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Escherichia coli*, *Salmonella enterica*, *Klebsiella pneumoniae*).<sup>1</sup>

Resistance arises due to mutations/insertions in genes involved in LPS biosynthesis (lpx, pmrA/B, mgrB, phoP/Q) and/or the acquisition of phosphoethanolamine transferases (PEtN). Of great concern is the recently described plasmid-encoded PEtN, MCR-1, now found worldwide in a range of animal, human and environmental bacterial isolates.<sup>1</sup>

CHROMagar™ COL-APSE is a sensitive and specific medium for the growth of Colistin resistant bacterial pathogens with a lower limit of detection of 10 CFU/mL. This new chromogenic medium may be useful as a primary isolation medium in the surveillance and recovery of Colistin resistant bacteria from complex human, veterinary and environmental samples especially those with plasmid mediated MCR-1 or novel mechanisms of polymyxin resistance.<sup>1</sup>

1- Novel Chromogenic Culture Media (CHROMagar™ COL-APSE) for the Isolation and Differentiation of Colistin Resistant Gram-negative Pathogens. Muhd Haziq F Abdul Momin, David W Wareham. ECCMID 2017.

### Medium Performance

#### 1 COLOURFUL DIFFERENTIATION OF COLONIES WITH ACQUIRED COLISTIN RESISTANCE

Distinguishing between *E. coli*, Coliforms, *Pseudomonas* and *Acinetobacter*.

#### 2 IMPRESSIVE LIMIT OF DETECTION

10 CFU/mL

### Medium Description

<b>Powder Base</b>	Total .....	42.5 g/L
	Agar .....	15.0
<b>+</b> <b>Supplement</b> (included in the pack)	Peptones .....	20.0
	Salt .....	5.0
	Chromogenic and selective mix .....	0.8
	Growth factors .....	1.7
	Storage at 15/30 °C - pH : 7.1 +/-0.2	
	Shelf Life .....	> 18 months
	Liquid form.....	2 mL/L
Storage at 15/30 °C		
Shelf Life .....	> 18 months	

Usual Samples	stools, rectal swabs, veterinary, environmental samples
Procedure	Direct Streaking. Incubation 18-24 h at 37 °C Aerobic conditions.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

### Order References


Please use these product references when contacting your local distributor:

5000 mL pack ..... CO262  
(included in this reference: powder base CO262(B) + supplement CO262(S))

Manufacturer: CHROMagar  
4 place du 18 juin 1940 75006 Paris - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)

Find your nearest distributor on  
[www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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