


For Research Use Only (RUO)
Not for use in diagnostic procedures

● CHROMagar™ O157




**For the selective isolation and differentiation
of *E.coli* O157**

CHR  **Magar**
The Chromogenic Media Pioneer

● CHROMagar™ O157

For Research Use Only (RUO). Not for use in diagnostic procedures.

For the selective isolation and differentiation of *E.coli* O157 in food samples




Medium Description

Powder Base	Total	29.2 g/L
	Agar	15.0
	Peptone and Yeast extract	13.0
	Chromogenic mix	1.2
	Storage at 15/30°C - pH: 6.9 ± 0.2	
Shelf Life	3 years	

Please read carefully the instructions for use (IFU document) available on www.drg-international.com

Plate Reading

- *E.coli* O157
→ mauve
- Other bacteria
→ steel blue, colourless or inhibited



Quality Control Strains

<i>E. coli</i> O157:H7 ATCC® 35150	mauve
<i>E. coli</i> ATCC® 25922	metallic blue
<i>Klebsiella</i> ATCC® 13883	metallic blue
<i>E. faecalis</i> ATCC® 29212	inhibited

ATCC® is a registered trademark of the American Type Culture Collection

Order References

Please use these product references when contacting DRG:

1000 ml pack EE220 5000 ml pack EE222 25 L pack EE223-25
Bulk on request

For US Customers: for Research Use Only

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France

Distributed in the US by DRG International
www.DRG-International.com