

# **DxDATA™**

Automated Chemiluminescent  
Immunoassay Analyzing System



**BioCheck**  
DRG Instruments GmbH

[www.DRG-Diagnostics.de](http://www.DRG-Diagnostics.de)

## FEATURES

**Accurate:** Magnetic particle chemiluminescent immunoassay method

**Convenient:** Single cartridge, low maintenance, no liquid path, easy to operate, easy to maintain

## SPECIFICATIONS

**Size:** 610 X 410 X 520 mm

**Weight:** 35 kg

**Sample throughput:** 8 channels

**Sample type:** Serum, Plasma, Whole blood

**Detection time:** 17-29 minutes

**LIS/HIS:** Unidirectional

## TEST MENU

### TBI (Traumatic Brain Injury):

GFAP  
PGP9.5

### Inflammation:

CRP  
PCT  
SAA  
IL-6

### Gonadal Hormone:

AMH  
Estradiol\*  
β-HCG  
FSH\*  
LH\*  
Progesterone\*  
Prolactin\*  
Testosterone

### Cardiovascular & Cerebrovascular Disease:

CKMB  
cTnI  
D-Dimer  
hsCRP  
Myo  
NT-proBNP

### Thyroid Function:

TSH  
T3\*  
fT3\*  
T4\*  
fT4\*  
Anti-TG  
Anti-TPO

### Diabetes:

HbA1c  
C-Peptide\*  
Insulin

### Bone Metabolism:

PTH  
25-OH Vitamin D